

Submission Form

Phone 03 489 4832 Email drl@drl.net.nz

DRL

Cre8ive9547 05/23

You can either print or complete on your computer. (remember to save the file once you have filled it out)

Include a print out of this form to send with your samples. Don't forget to email drl@drl.co.nz

to book your samples in, as priority will be given to booked ones.

Customer Details		
Farm		Vet
Farm Manager:		Vet Name:
Farm Name:		Vet Clinic:
Email (for invoice):		Email (for invoice):
Address:		Address:
Sample		
Species: Cattle Deer Sheep Goat Other		
Sample: Blood Faeces Tissue Milk Number of Samples:		
Tests required		
Serum	Faecal PCR	Ear Notch Milk
Johne's Disease	Johne's Disease (Individual)) BVD Johne's Disease
ELISA + Paralisa™	Johne's Disease (Pooled)	BVD Calf < 35 days BVD
BVD		CAE
CAE		
Raw Data (surcharge of \$0.50 +GST per sample)		Date Required:
Comments:		

Please label all samples clearly and courier samples overnight to:

DRL ^c/ Invermay Research Centre, 176 Puddle Alley, Mosgiel 9092

Samples booked in will be afforded priority. Please email this form to drl@drl.net.nz

www.drl.net.nz

Notes and Pricing

Please contact the lab if you need assistance sourcing any sample supplies 03 489 4832 or drl@drl.net.nz

Serum Samples:

- Blood (red or purple top tubes; min. 3ml).
- While refrigeration is desirable it is not essential; please protect samples from extremes of heat or cold.
- Do not freeze.

Milk Samples:

- Please use screwtop plastic containers for milk samples. Min. 3ml per sample.
- Please keep refrigerated and pack with ice packs in a chilly bin.

Individual Faecal PCR: \$40+GST*

- Please use screwtop plastic containers for dung.
- Use a fresh glove to collect each sample and avoid cross contamination of samples.
- Min. 3g per sample, ideally 5-10g (1 tablespoon approx.)

Pooled Faecal PCR: \$50+GST* (Pools of 10)

- If sample pooling is requested, send individual samples as we will carry out the pooling here, that way we can easily retest individual samples from positive pools as required.
- If a pool is positive, we will contact you concerning your preference for individual sample re-testing.
- Please use plastic screwtop containers for dung. Use a fresh glove to collect each sample and avoid any cross contamination of samples.
- Use a fresh glove to collect each sample and avoid cross contamination of samples.
- Dung samples freshly deposited on the ground are acceptable as long they can be unambiguously linked to their originator and land on an otherwise clean surface.

*special pricing is available for custom diagnostic programmes, please call the lab 03 489 4832 to discuss or Simon on 021 249 7710.

Terms and Conditions

DRL retains ownership of all information and intellectual property derived from the analysis of samples submitted.

All invoices are to be paid by the 20th of the month following (normal commercial practice).

If any debt remains unpaid after 90 days, it will be referred to a collection agency. Collection costs will be incurred by the debtor.

